HELP FOR THE NEW SPECTROMETER UNDER DELPHI 4.0 (2001)

Julien Bobroff, Safia Ouazi, Véronique Brouet
LPS, Bat.510, Université Paris XI, 91405 Orsay Cedex, France
bobroff@lps.u-psud.fr

How to start the program.....	2
How to work with the different windows :	2
Where is the data saved and how does it look like ?	2
Some bugs you may (will) meet :.....	3
The acquisition RMN window.....	4
The TFRMN window	6
The Pulse Sequence window	8
Structure of the file c:\newspectro\Parametres RMN.txt	12
interpréteur window:.....	13
Syntax for programming	13
The variables which cannot be used for everything:	14
The instructions	15
Instructions for acquisition parameters	15
Instructions for the boite noire	16
Instructions for analysis of the data	16
Instructions for manipulation of an array	16
Some examples of programs	17
Reception Amplifiers Window.....	20
Variation of the amplification versus the different gain values :.....	20
Help for those who would like to program	23
General rules :	23
Data Format	23
Structure of the program	23
VERONIQUE's NOTES.....	28
Pannes de spectro	28
Initial setup when installing the program	34
Initial setup when installing the program	34

How to start the program

start Delphi Borland 4.0 (the yellow greek temple icon). Then press Ctrl+F11 and choose in the directory c:\newspectro\pilotageManip3 the project named PilotageManip. Then press F9 (being sure that the spectrometer is open).

How to work with the different windows :

All the acquisition is done through the **acquisition RMN** window. The measured data (sum of all the shots) is in ECHO which appears in A/B of this window. The last single shot data is in C/D in OSCILLO. One uses the Edit box of this window (in A/B) to save the data or recall another one (in C/D) for comparison purpose.

There is an automatic soustraction of background (taken between 75% and 100% of your data) on ECHO at the end of the run or at each rafraichissement. (no soustraction like that on OSCILLO).

All the Fourier Transforms are done in **TFRMN** window. Everything works as in the acquisition window, with the TF parameters appearing in front pannel. The TF can be done on data acquired in the acquisition window, or on previously saved data.

There is an **interpreteur** window aimed at more sophisticated commands or list of commands. It can be used for automatic runs, or for automatic analysis of the data. The syntax to be used is not easy, so it is easier to use the preexisting examples.

There is a **reception** window where to find and change the spectrometer reception amplifier parameters.

There is a **sequence** window which contains the pulse sequence. Do not be confused with the interpreteur window : the sequence window is specifically dealing with the radiofrequency pulses sent to the sample. It has a list of parameters on its left, and a box where to write the sequence on its right. The sequence uses the parameters of left part.

There is a **gestgage** window wich contains the acquisition card (GAGE card) parameters, like digitilization, voltage scale, number of points, etc.

Where is the data saved and how does it look like ?

All data are saved in the current directory, which appears in the acquisition RMN window. Typical saved data have the form : blablaECHO1.spe, blablaECHO2.spe, blablaTF3.spe, blablaECHO4.spe, etc... The number of the last saved data appears near the directory name. You can also save data in specific names as bloubli.spe. These .spe files are the program format files and can be read only by this program. You can also save your data in additional files in ascii format to export them elsewhere (in .dat extension).

An important constraint : at beginning of the program, this directory is found by the software in an ascii file named current.txt located in c:\newspectro\data\ containing just the full directory name (like c:\newspectro\data\blabla). The software also needs a file named

info.txt located in the current directory (here c:\newspectro\data\blabla) which contains in the first line again the directory name (here c:\newspectro\data\blabla) and in the second line the number of the last saved file (here, it would be 4). If one of these files does not exist, you have to create it.

What to do if a bug occurs and stops the program ?

Try to switch to the Delphi window (the yellow greek temple icon) with Alt+Tab. Then, press the F9 key to make it run again. It may bug again, say always OK and F9 again. Switch back to the program window with Alt+Tab when it works. If it always fails and does not start again, then press Ctrl+F2 to stop the program. Start it again with F9.

Some bugs you may (will) meet :

When starting, the program does not find info.txt or current.txt files or c:\newspectro\Parametres RMN.txt (see paragraph Where is the data saved) : it will crush.

The rf pulses or frequencies don't seem to work : try to change the value of the working frequency in the reception window (at 14000000) to 14000001 and put it back to 14000000.


When starting, there is a division by zero mistake : no idea where it comes from, doesn't matter.

Problems when trying to call back a spectrum : check in the current directory if it really exists

Problems with the interpreter commands : did you really use brackets, "" for characters, and ';' at the end of the line (like command["aa",1,2];) ?

Problems with the sequence window : did you really put a dot '.' at the end of the program ? Do you really put time in msec ?

The acquisition RMN window

in the Edit box, you can type the following commands :

s : standard saving of the data. It saves the data currently in A/B channel in the program internal format (which is then readable again by this program but no other). Each data is saved at c:\NewSpectro\data\example\exampleECHO4.spe if "example" is the current directory, and if already 3 files were saved before (either TF or ECHO). In the box of data A, only exampleECHO4 will appear.

sblabla : saving with a specific name (here blabla). It saves the data in A/B at c:\NewSpectro\data\example\blabla.spe (instead of the preformatted name exampleECHO4)

ascii : saving of the data in ascii format. It saves A/B in an ASCII format readable with other programs like Exel or Origin. The data consists of 3 columns : x, y of channel A, y of channel B. If the present name in the box appears for example as exampleECHO4, it will save it as c:\NewSpectro\data\example\exampleECHO4.dat
Caution : the ascii data cannot be recalled then in this program and is usually bigger in size than the .spe files.

asciblaba : saving of the data in ascii format with a specific name (here blaba). It acts like sasc but saves it as c:\NewSpectro\data\example\blaba.dat

r4 : recalls data (here n°4) which was saved previously with s and puts it in C and D.

rblabla : recalls data with specific name (here blabla) which was saved previously with s and puts it in C and D.

add5,7,8 or add5-8 : adds the files 5,7 and 8 or 5 to 8 in SOMME

bg or bg3 or bgnom : subtracts background on A/B or 3 or nom (taken between 3/4 and end of the data)


int(A,5,24.3) : integral of channel A between 5 and 24.3

intn(A,5,24.3) : integral of channel A between 5 and 24.3 normalized = average value per channel

intecho(A,5,24.3,50,100.3) : echo intensity = square root of [square of (integral of channel A after subtraction of bgd taken between 50 and 100.3) + square of (idem for B)] . Note that you don't need to take the same interval width for background and echo max. If you put C, it will take C and D. If the two background values are equal, it does not perform the bgd soustraction.

intechon(A,5,24.3,50,100.3) : echo intensity normalized = average value of the echo per channel

The TFRMN window

in the Edit box, you can type the following commands :

Like in Acquisition RMN window, you can use **s**, **sasc**, **r**. Here, the data will be saved as exampleTF4.spe instead of exampleECHO4.spe.

tf : does the Fourier Transform of A/B of Acquisition RMN window. It uses the parameters of Fourier Transform displayed in the window and puts the result in A/B with name TF.

tfblabla : does the Fourier Transform of the file blabla.spe.

tf5 : does the Fourier Transform of the file with number 5 (in our example a file of the form c:\NewSpectro\data\example\exampleECHO5.dat)

tf5-9 : does the Fourier Transform of the file with number 5 to number 9 and saves them each time (with same parameters for all TF)

The Fourier Transform Parameters : they can be saved and recalled by the button. They consist of ascii file with name like blabla.par and with explicit form (you can read them with word for exemple).

Time Zero : the origin of time of the signal in microseconds

Binning : it averages the time data by replacing a group of binning points by its average. Put 1 doesn't do any binning. You can then put 2, 3, etc... to see the effect. It smoothes the time signal. If you bin too much, you loose some informations contained in the signal.

BGD : the interval (in microseconds) where to take the background which will be subtracted

Echo Length : the length of the echo in microseconds to be used (if zero, it uses all the data)

Nb Points for computation : you can use more than the time signal you have to do a more smooth TF by adding some fake points with zero intensity after your signal (Filling Zero). You do so by choosing in Nb of Points for computation a higher number of points than the ones you really recorded. The maximum is a priori the better (16384).

Nb Points to keep for TF : usually, you are interested only in the center of the TF, not in all the points (like the 16384). You can then restrict your TF to this number of points (the program just destroys the other points), it creates a much lighter file easier to work with. Use typically 1000.

Frequency Step : the frequency step between two TF points. It is equal to $1/(\text{Nb of Points for computation} \times \text{Time step} \times \text{binning})$ and is given in kHz.

Example of a time signal and the parameters to be chosen, and the corresponding TF :

The Pulse Sequence window

Click this menu to save/load sequence

Name of sequence

Parameters used the sequence (time in msec, phase in 0.1°)

Program of pulse sequence (here, $\pi/2$ -tau- π and $\pi/2$ inverted-tau- π , with turning reception)

Click to change the acquisition mode (here, the reception is turning, so no need for clicking)

Status of the pulser

Click to start/stop the pulser (it will automatically start if you start an acquisition)

you can recall an old sequence. But take care : if you modify it and then start the sequence, it saves the new version. So keep somewhere a backup with another name of your standard sequences. You can enter by hand the parameters on left side. It consists of 3 columns :

Phase variation : if standard echo with all phases at 0° gives A,B, then when varying the phases of the different pulses, here is the result :


Phase $\pi/2$	Phase π	Phase Reception	Channel 1	Channel 2
0°	0°	0°	+A	+B
180°	0°	0°	-A	-B
0°	0°	90°	+B	-A
0°	0°	180°	-A	-B
0°	0°	270°	-B	+A

A simple $\pi/2$ - tau - π :

$\pi/2$ length = TA with gain AA, 0°
 Tau = TC + $3\ \mu\text{s}$ (due to blanking)
 Pi length = TB with gain AA, 0°
 Rep time = TR (in fact, = TR + 3 + TA + TC + TB + TC + 2)
 Blanking on Channel Out E
 Trigger for Gage on Channel F
 Trigger for scope on Channel Out A
 Pulses on channel Out Pulse

T50,CA; starts a trigger (for scope use) on CA during $5\ \mu\text{s}$
T30,CE,PB01,GA0; starts a blanking pulse on CE during $3\ \mu\text{s}$
TA,CE,GAAA; starts an rf pulse with amplitude AA during time TA (and keeps on pulse CE)
TC,GA0; waits during TC with no pulse
T30,CE,PB01,GA0; starts a blanking pulse on CE during $3\ \mu\text{s}$
TB,CE,GAAA; starts an rf pulse with amplitude AA during time TB (and keeps on pulse CE)
T20,CF,GA0; gives a $2\ \mu\text{s}$ trigger on channel F
TC,GA0,PB01; waits during TC with no pulse
TR,PB01,GA0; waits during TR with no pulse
J. end of sequence

hereis the corresponding histogram:

A $\pi/2$ - tau - π alternated with $\pi/2$ inverted- tau - π :

Two ways : either make reception turn as well in phase to have the recombination of pulses, or make recombination in the program.

Advantage and inconvenients : if recombination by turning reception, allows a fast repetition but does not cancel well the saturation after π pulse (because this saturation has not been turned in phase). If recombination by soft, obliges a slow repetition not to miss any shot, but cancels well the saturation.

To activate the program mode recombination (if your reception phase is left 0°), you have to click the mode in the sequence window. It will do it automatically if in the name of the sequence, you put somewhere 'mode 2' or 'mode 4'.

For example, a sequence named 'piSur2 pi mode 2.seq' when loaded will automatically activate the mode 2 in the window.

If you forget to activate the mode when needed, your resulting data will always be zero !

2 Steps Method with recombination by turning reception :

```
T30,CF;
T30,CE,PB01,GA0;
TA,CE,GAAA;
T50,CF,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB01;
TR,GA0,PB01;

T30,CF;
T30,CE,PB1801,GA0; here pulse pi/2 is inverted
TA,CE,GAAA;
T50,CF,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB1801; here reception is turned
TR,GA0,PB1801;
J.
```

2 Steps Method with recombination with reception always at 0°:

One must click then on Mode Acquisition 0° 2 pulses option and choose a slow repetition time

```
T30,CF;
T30,CE,PB01,GA0;
TA,CE,GAAA;
T50,CF,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB01;
TR,GA0,PB01;

T30,CF;
T30,CE,PB1801,GA0; pi/2 inverted
TA,CE,GAAA;
T50,CF,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB01; reception is not turned
TR,GA0,PB01;
J.
```

4 Steps Method with recombination by turning reception :

```
T30,CF;
T30,CE,PB01,GA0;
TA,CE,GAAA;
T50,CF,GA0,PB01;
```

TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0;
TR,GA0,PB01;

T30,CE,PB901,GA0;
TA,CE,GAAA;
T50,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB2701;
TR,GA0;

T30,CE,PB1801,GA0;
TA,CE,GAAA;
T50,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB1801;
TR,GA0;

T30,CE,PB2701,GA0;
TA,CE,GAAA;
T50,GA0,PB01;
TC,GA0;
T30,CE,PB901,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB01;
TR,GA0;
J.

4 Steps Method with recombination with reception always at 0°:
One must click then on Mode Acquisition 0° 4 pulses option and choose a slow repetition time

T30,CF;
T30,CE,PB01,GA0;
TA,CE,GAAA;
T50,CF,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0;
TR,GA0,PB01;

T30,CE,PB901,GA0;
TA,CE,GAAA;
T50,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;
TB,CE,GAAA;
T50,CA,GA0,PB01;
TR,GA0;

T30,CE,PB1801,GA0;
TA,CE,GAAA;
T50,GA0,PB01;
TC,GA0;
T30,CE,PB01,GA0;

TB,CE,GAAA;
 T50,CA,GA0,PB01;
 TR,GA0;

 T30,CE,PB2701,GA0;
 TA,CE,GAAA;
 T50,GA0,PB01;
 TC,GA0;
 T30,CE,PB01,GA0;
 TB,CE,GAAA;
 T50,CA,GA0,PB01;
 TR,GA0;
 J.

Structure of the file c:\newspectro\Parametres RMN.txt

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1.00

DC

Last Pulse Sequence (exple : nom.seq (found in sequenceRMN))


C:\NewSpectro\sequenceRMNecho 2 pulses mode vary.seq

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Status of the program

If you want to execute a single line command, put it here and return

To access this menu (for load/save the script), press right button of the mouse

The script of the program (here a pi/2-pi/2 search with plot of the result)

It allows to program automatically a lot of instructions, either for doing automatic experiments or also automatic analysis of the data.

Syntax for programming

The syntax for programs is very simple :

- each line must end with ';'.
- you can use variables. They must be in capital letters. For example :

RESULTAT=5;

This will affect the value 5 to the variable RESULTAT. All variables are real numbers (no string no boolean). There are some specific names of variables you cannot use for everything (which are listed below). You can follow the value and name of all variables by clicking on the menu 'variables' when you are in the interpreteur window.

- you can use labels following the syntax :
instruction 1;
instruction 2;
la: instruction 3;
instruction 4;
goto la;

here, instruction 3 and 4 will be in an endless loop. Note that there must be an instruction on the line 'la:...'. You can use for labels the name la, lb, lc, ld, le ...

- you can use conditions if - then following for example:

```

NUMBERRUN=1;
La:instruction 1;
Instruction 2;
NUMBERRUN=NUMBERRUN+1;
If NUMBERRUN<= 10 then la;
Instruction 3;

```

Here, the program will do instruction 1 and 2 ten times, and then go to instruction 3

- you can add comments after ';' like (very usefull next time you look at your program)
- you can use arithmetics functions (all are listed in the menu 'fonctions' when you are in the interpreteur window) for example :

```

X=5;
RESULTAT=ABS(-COS(X)); this will put in resultat absolute value of -cos(5)

```

The variables which cannot be used for everything:

Some names of variables are booked for specific tasks. For example, when you write 'TEMPERATURE=200.5;', It will not just put 200.5 in a variable TEMPERATURE, but It will also program the real temperature through ITC4 to go to 200.5. So take care not to use the following names for any purpose, you would have big big surprises !!!

Variable (and an exemple of instruction)	What it does
NB=5000;	Puts the number of points to do for acquisition to 5000
TEMPERATURE=300.2;	Puts the temperature in ITC4 to 300.2, putting at the same time heater and gaz flow automatic
TA=0.002;	Puts the pulse TA to 0.002 msec (idem for TB,TC,TD,...)
AA=2000;	Put the amplitude AA of pulse to 2000 (idem for AB,AC,AD...)
PA=1801;	Put the phase of pulse PA to 180° (idem for PB,PC,...)
TZERO=500.54; BIN=3; BGDINF=1254.3; BGDSUP=1354.6; LEFTLENGTH=150.65; RIGHTLENGTH=150.65; NBCOMPUTE=4096; NBTF=1000; FREQ=43252.25	Changes all these parameters for Fourier TTransform

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Note that there are already many instructions listed above in the array for specific variables.

All instructions have the form :

Instruction; here, instruction with no input or output variable

Instruction[24,35.5]; here, instruction with two real number variables

Instruction["RESULTAT","ZZ",24.5]; here, instruction with two string variables and one real number variable.

Here are some instructions :

Instructions for acquisition parameters

Startnosave; starts an experiment but does not save it at the end. It will appear as ECHO in the acquisition window.

Start[1]; starts and saves at the end the run in the usual syntax

Start[5]; starts and saves at the end the run in the usual syntax five times. Note that the user wants to interrupt this instruction, if he presses stop from the acquisition window, it will stop the run but do all the others. You need to interrupt delphi itself !

Start["blabla"]; : idem but saves the run under the name blabla

Save["blabla"]; saves the data named blabla in c:\newspectro\data\Echantillon\Echantillonblabla5.spe.

Save["blabla","blabBis2"]; saves the data named blabla (for exple, ECHO at the end of a run) in c:\newspectro\data\Echantillon\blablaBis.spe.

saveascii["nom","nom2"]; saves nom in ASCII 3 columns in nom2.dat in the working directory

ChangePulse[« TA » ,0.01] ; changes a pulse value. Here, it changes TA to the value 0.01 in SeqPar window. Caution : it does not mean that the pulse has been physically changed; you have to type then startSequence to validate this change.

CompteNbCoupsDejaFait["NBCOUPSDEJAFIT"]; puts in the variable NBCOUPSDEJAFIT the actual nb of shots already done, not the one to be done (usefull if the user stops the experiment from the acquisition window during the unfinished run to interrupt before NB has been done)

AddAndReplace["ECHO","NOM"]; : it will recall NOM.spe, will add to it ECHO and will save back NOM.spe. If NOM.spe does not exist, it will create it. If ECHO does not exist, it will do nothing.

AddAndReplaceT1["T1300K","TEMPS"]; if for exple the value in TEMPS is 0.5, and the name of the working directory is c:\NewSpectro\data\direc, it will recall direcT1300K0_5.spe, and will add to it the last data in ECHO, then saves it back. If direcT1300K0_5.spe does not exist, it will create it

Instructions for the boite noire

AddComment["mesure en cours"] : adds the comment mesure en cours to the black box
AddCommentVariable["ZZ"]; adds the comment 'ZZ=50' if the variable ZZ exists (50 being its value here)

AddCommentVariable["ZZ","commentaire quelconque"]; idem and adds at the end of the line 'commentaire quelconque' which is a free comment

ParamGAGE : puts in the boite noire the parameters used by GAGE card

ParamReception : puts in the boite noire the parameters used for reception (amplis)

Instructions for analysis of the data

Integ["blabla",1,15.5,25.2] : integrates channel 1 of data blabla between 15.5 and 25.2. Result is zero if there was a mistake.

TF[5] : Fourier Transform of file 5 using the parameters of the TF window

TF["blabla"] : Fourier Transform of file blabla.spe using the parameters of the TF window

TF["blabla","parametre"] or **TF[12,"parametre"]** : Fourier Transform of file blabla.spe or file n°12 using the parameters of the contained in parametre.tfp saved previously in the TF window

integEchoT1[10.3,14.3,25.2,26.2,"INTEGRALE"] : puts in INTEGRALE the value of the integrale of the two channels of last data saved where signal is taken between 10.3 and 14.3, after substraction of background between 25.2 and 26.2.

integEchoT1["LEFT","RIGHT","BGLEFT","BGRIGHT","INTEGRALE"]; exactly the same but takes the value of the variables LEFT, RIGHT,... for the min and max of cursors.

Instructions for manipulation of an array

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Some examples of programs

A simple experiment

Addcomment["experiment at 300K"]; adds the comment in the black box
 NB=20000; specifies nb of points for each acquisition
 TR=50; repetition time in the sequence (msec)
 Start[1]; starts and saves one acquisition
 NB=30000; specifies another nb of points for each acquisition
 Start[5]; starts and saves five acquisitions

A search for $\pi/2$ -pi optimal

Addcomment["loop for $\pi/2$ -pi search"]; adds the comment in the black box
 NB=20000; specifies nb of points for each acquisition
 TR=50; repetition time in the sequence (msec)
 PulseIncrement=0.001; increment in msec between two steps for the search
 TADebut=0.001; starting pulse length for $\pi/2$
 TAFin=0.01; ending pulse
 TAFin=TAFin+PulseIncrement; in order to do the last step also
 TA=TADebut; initialize the value for $\pi/2$
 TB=2*TA; idem for pi
 la:start[1]; begining of the loop : starts and saves the acquisition
 TA=TA+PulseIncrement; at each step, we increase TA
 TB=2*TA; we increase TB
 if TA<TAFin then la; end of the loop : if we reach the last length, we stop the loop
 TA=TADebut; reinitialization of the pulse length for $\pi/2$
 TB=2*TADebut; reinitialization of the pulse length for pi

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saturate for an input signal of about -35 dB, the signal is deformed (but it is not easily detectable).

GainLo0,GainLo90 : Gain applied to the Lo frequency. The total gain should not depend from these values as long as the RF level is sufficient for the demodulation to work properly. In practice, values between 100 and 200 are correct and 181 is usually used. There are sometimes problems with offsets associated with these gains (see below *how to set gains*).

GainBF0,GainBF90 : Two low frequency gain (40 dB) on separate channels. For them, minimum corresponds to 0 and maximum to 255 (see graph).

How to set gains:

The usual values are :

GAINHF1 = 140	GAINHF2 = 160	At 75 MHz : gain between 60 and 70 dB
GAINLO0 = 181	GAINLO90 = 181	noise : 200 mV with 190
GAINBF0 = 170 to 190	GAINBF90 = 170 to 190	

One cannot set gains independently from each other. For given values of (GAINHF1,GAINHF2) (chosen to avoid saturation), there is only a limited range of (GAINBF0,GAINBF90) and (GAINLO0,GAINLO90) values that will work properly, which means they don't deform the signal, keep the quadrature to 90° and correspond to identical gains on the two channels. I usually keep GAINHF1 and GAINHF2 fixed and change GAINBF0 and GAINBF90 from 170 to 190 depending on the signal amplitude.

The choice of (GAINLO0,GAINLO90) is crucial to minimize **offsets**. With a wrong setting, these offsets can be large and dissymmetric on the two channels (which can result in slightly different gains for the two channels even for the same parameter values). They are sometimes unstable with time, which can cause baseline problems (especially at low frequency ?). These offsets are not subtracted with the usual sequences, see *acquisition* to see how to get rid of them and/or compensate unequal gains on 0° and 90°, if necessary. For some frequencies, these offsets obviously oscillate (look for example at 70 MHz). This occurs for multiples of 14 MHz, so that it might be necessary to shift the reference 14 MHz in order to work at these frequencies (but care must be taken that the operating frequency for the reception (normally 210 MHz =15*14 MHz) is still OK).

Other settings :


- For a higher gain (about 80dB which is near the maximum):

GAINHF1 = 140	GAINHF2 = 120
GAINLO0 = 131	GAINLO90 = 131
GAINBF0 = 200	GAINBF90 = 200

- For low frequency (used by Feri)


GAINHF1 = 175	GAINHF2 = 161
GAINLO0 = 208	GAINLO90 = 204
GAINBF0 = 184	GAINBF90 = 182

Warning : the gain seems to be anomalously low below 20-30 MHz. To be checked.


This plot shows the variation of the gain of the reception when changing hf1,hf2

(Signal taken at 75 MHz, with a signal -60dB, with the values
GainHf1/Hf2=140/160;GainLo 0/Lo90=181/181;GainBF0,BF 90=170/170;)

This plot shows the variation of the gain of the reception when changing Lo0,Lo90

(Signal taken at 75 MHz, with a signal -60dB, with the values
GainHf1/Hf2=140/160;GainLo 0/Lo90=181/181;GainBF0,BF 90=170/170;)

This plot shows the variation of the gain of the reception when changing BF0,BF90,

(Signal taken at 75 MHz, with a signal -60dB, with the values
GainHf1/Hf2=140/160;GainLo 0/Lo90=181/181;GainBF0,BF 90=170/170;)

Help for those who would like to program

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General rules :

when you want to change something, first do a backup on the disk of the present version of the program, and put the date of the backup in the title of the directory (exple : c:\newspectro\backup pilotagemanip3 11 august 2001).

When you add something, put a remark as {groupe Alloul } in front of any new procedure. Try to put it also at the end of the unit if possible.

Data Format

The data format and its use is done in **editdata** :

All data are contained in a variable named und.spect which has the following properties :

Und.spect[i].nom is the name (exple 'ECHO') of data number i

Und.spect[i].X is the number of points

Und.spect[i].Y must be put at 1 (one dimensional data, no Y scale)

Und.spect[i].comment is a possible comment (exple 'this is a signal of NMR')

Und.spect[i].typex is defining the type of x data to distinguish between different kinds (here 'rmn')

Und.spect[i].typey is here single (no 2nd dimension)

Und.spect[i].data[j,0].x contains the x value of the jth point of the data (second number at 0 no as 2nd dimension, so that Und.spect[i].data[j,0].y is not used)

Und.spect[i].data[j,0].z contains the z value

Und.spect[i].data[j,0].mon contains another possible value

For NMR time signal: we use .x for the time, .z and .mon for real and imaginary intensity

For NMR Fourier Transform: we use .x for the frequency, .z and .mon for real and imaginary intensity

Indexof(name) gives the index associated to name. For example, to use the data 'ECHO', you can use und.spect[indexof('ECHO')]

One can also use different procedures applicable to und.spect :

Und.spect.clear, total, minmax, savetofile, readfromfile, savetoPAXY, readfromPAXY, free, procedure free; addspectre, deletespectre, addition, soustraction, saveall, add, addpt.

Structure of the program

Important units/Files are the following:

* **PilotageManip** : the source of the project, which contains which window is created at beginning
(exple : Application.CreateForm(TseqPar, SeqPar);) and which units are used in the full program.

* **PilotageManipMainForm/PrincipalForm** :the main window which contains all the other. Contains the initialisation of the initial data :

Exple :

```
ind:=und.addspectre('DATA',4096,1);
und.spect[ind].comment:='RMN';
und.spect[ind].typex:='rmn';
und.spect[ind].typey:='single';
for i:=0 to 4095 do begin
  und.spect[ind].data[i,0].x:=i;
  und.spect[ind].data[i,0].z:=1000*sin(i*2*3.14/10)*exp(-i/1000);
  und.spect[ind].data[i,0].mon:=1000*cos(i*2*3.14/10)*exp(-i/1000);
end;
```

At beginning, it creates 4 files : DATA, ECHO, OSCILLO, TF of the rmn type.

* **AcquisitionRMN/NMRForm** : a graphic window which is aimed at the acquisition of the data only. Thus, the appearance of data is a little bit different than in usual graphic window : each data appears in two channel at once (A,B or C,D), and the .z part appears in A, the .mon part appears in B (respectively C,D). So real and imaginary part of the signal are always appearing together.

There is a Directory window and a n° window. These two contain the working directory and the number of the last saved data. To make it work, the program uses

- a file named c:\newspectro\data\current.txt which contains in text format the full path of the working directory (exple : c:\newspectro\data\nouveau). If this file is not present, the program cannot work.
- A file in the working directory named info.txt which contains again the working directory and the number of last saved filed

exple :

c:\newspectro\data\nouveau

56

If this file is not present, again it will not work.

To change working directory, press the directory button. Note that these informations are the same in the TF window as well.

The code of **acquisitionRMN** is modified as compared to a standard graphic window : it contains now a modified version of **tfrspectre** to allow the special im/real part display explained hereabove.

It is also linked to real acquisition of data through **Rafrachir** : a procedure called during acquisition and linked to a timer which regularly refresh the A/B channel putting in it the total signal, and in C/D the last instantaneous signal. So A/B is called ECHO, and C/D OSCILLO.

It has a blank command box above the graph which allows specific NMR commands. These commands appear in **traitchaine** in the following way: if for exemple to save a file, you want to type s or S you write the code in traitchain :
if (s='S') or (s='s') then

```

begin
  dataeditform.sauveNMR(cocourbe1.text); {this goes in editdata and uses the procedure
sauveNMR to save data}
  RafrachissementPileEtDirectoryNMR; {this calls the refreshment of number of the
spectrum and directory informations}
  cocourbe1.text:=und.spect[1].nom; {this puts the new data name in the front pannel}
end;

```

You can add there your specific commands for this window.

* **CompilRMN** : contains especially **compilo** , the procedure which compiles a pulse sequence from the SeqPar window. This is where the syntax of compilation of sequences is written.

* **EditData** : as explained above, contains the procedures and functions which deal with the data contained in und.spect. We added there specific procedures for NMR data format :

- sauveNMR saves the data in A/B. If the data is named ECHO, and the current directory is c:\newspectro\data\nouveau and the last saved file has number 10, it will save it in the name nouveauECHO11.spe in a format readable only by the program. This is similar with any other name (TF will lead to nouveauTF11).
- recallNMR recalls the data saved through sauveNMR (i.e in the format nouveauECHO10.spe if directory is named nouveau)
- DecaleSpectres : used in sauve and recall, reorganises the spectra in spectlist :
1st data is the present data, either the last one recalled or saved, or acquired...
2nd data is ECHO
3rd data is OSCILLO
4th data is TF
5th data is the previously called or saved,
6th data is the previously called or saved,
etc... up to 10 data maximum (this limit number of authorized files in spectlist is fixed in decalespectres in the variable MaxiNombreDeFichierDansSpectListe). - -
- TransfertSpectre : used by DecaleSpectres, please do not use it from anywhere else, it is very tricky.

* **Interpreteur/InterpreteurForm** : the unit allowing batchs for programation. If you want to add your own command, do it using one already existing as a model. For exemple, here is the function which changes a parameter in the pulse sequence window seqPar : one uses for exemple ChangePulse["Ta",0.01];

```

function ChangePulse:integer;far;
begin
  SeqPar.ChangeParametreSequence(pars[1],par[2]);
  result:=1;
end;

```

one must also add to the list of commands ordre.addobject the line :

```

ordre.addobject('ChangePulse',tobject(@ChangePulse));
Each command has a form like : name["blabla",0.03,0.04,"boubou"];

```

Here, name is the name of the function called, and each variable is either a string (here blabla and boubou) or a real number (here 0.03 and 0.04). When you use these variables,

the string are noted **pars** and the numbers **par**. So here `pars[1] = 'blabla'`, `par[2] = 0.03`, `par[3] = 0.04`, `pars[4] = 'boubou'`.

If you want to add a variable that you use then (example the variable ALPHA), first go to the line

```
if NMRFORM<>nil then begin
  addvariable('TA',0,0,0,0,"",tobject(@ChangePulse));
  and put the line
  addvariable('ALPHA',0,0,0,0,"",tobject(@rien));
```

A more sophisticated command for variable allows you to use in the interpreter a line like ALPHA=4. Let's assume we want that this command makes the alpha box in the TFNMRForm Window to be then changed to 4. The way to do it is : do the addvariable as explained above, but with the line

```
addvariable('ALPHA',0,0,0,0,"",tobject(@ChangeAlphaValue));
```

Then, add a function as explained above of the form :

```
function ChangeAlphaValue(s:string):integer;far;
var Sender:TObject;x:extended;chaine:string;
begin
  x:=param[variable.indexof(s)].value;
  str(x:10:2,chaine);
  TFNMRForm.AlphaBox.text:=chaine;
  result:=1;
end;
```

Finally, do not forget to add the command at the proper place

```
ordre.addobject('ChangeAlphaValue',tobject(@ChangeAlphaValue));
```

- **GestGage/acquisParam** : the unit governing the gage card for acquisition of the NMR signal.
- **ParSeq/SeqPar** : the unit governing the pulse sequences.
- **Reception/Fenetre_Reception** : the unit governing the reception parameters (amplifiers and frequency of the reception part of the spectrometer).

TFRMN/TFNMRForm : exactly like **acquisitionRMN**, but for TF purpose. All the S, R, directory procedures are the same. But it is also designed to do TF

When you add a new window/unit : It must be open through the pilotagemanip unit as

```
Application.CreateForm(TSeqPar, SeqPar);
```

You have in the Inspecteur d'objet of the window to put on the parameter FormStyle the choice fsMDIChild. You also have to put a TAG number not used by the other windows (each window must have its own TAG). TAG allows then the database to remember the parameters of the window.

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